## Rec-P25 (MVV)

Vial of 100µg / Flacon de 100µg

Ref. VP002A

**STORE AT**  
**2-8 °C**

For research use only

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**CONSERVER A**  
**2-8 °C**

**Origin:**

Rec-P25 (MVV), purified and stabilised. Rec-P25 is a recombinant gag protein derived from Maedi Visna Virus (MVV), specifically targeted by antibodies which may be developed in infected sheep or goat.

**Presentation:**

Vial containing approximately 100 µg of rec-P25 lyophilised in presence of additives and preservatives.

**Reconstitution:**

Each vial must be restored with 0.2 ml distilled water in order to obtain a preparation at 500µg/ ml. It can then be diluted using the appropriate solution.

**Excipients:**

– Phosphate, Mannitol, NaCl.

**Purity:**

Purified recombinant P25, with a major band of about 26 ± 4 KDaltons on SDS-PAGE. Presence of P25 dimers can occur.

**Safety**

Any product of biological origin must be handled with all the required cautions, as being potentially infectious.

**Use:**

For in vitro use only.  
All research studies and protocols where a source of purified rec-P25 is necessary.

**Stability:**

**Restored :**

- 24 hours at room temperature
- 72 hours at 2-8 °C
- 2 months at -20 °C or below

**Lyophilized :** - Until the expiration date printed on the vial.

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CONSERVER A  
2-8 °C

## ANALYSIS CERTIFICATE

ANALYTICAL DATA	SPECIFICATIONS
1. <b>Protein Content</b>  ..... µg/vial	> 90 µg/vial
2. <b>SDS-PAGE (12 % acrylamide)</b>  1 major band of about ..... Kdaltons	1 major band of about 26 ± 4 Kdaltons

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DATE :

CONCLUSIONS :  Passed/Conforme  Refused/Non-conforme

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